



Endoscopic Bedside Kits

Removing visible bioburden after a procedure is a critical first step in the reprocessing of a flexible endoscope. The First Step Kit has simplified this process by providing the supplies needed to begin pre-cleaning a flexible endoscope immediately following a procedure. All Kits contain an endoscopic cleaning pad for wiping down the insertion tube and Simple2™ Enzymatic Detergent contained in a convenient pouch. Everything needed for immediate pre-cleaning before the scope is transported for reprocessing. First Step Bedside Kits are available as Ready-to-Use Kits or Add Water Kits and are available in sizes: 100 ml, 200 ml and 500 ml. See page 2 for details.



First Step™

- 100% compliant with SGNA and manufacturer guidelines.
- Safe for use with all flexible endoscopes.
- Offers a variety of cleaning pad styles.
- Ensures the proper dilution of enzymatic detergent.

The Ready-to-Use First Step Kit

Includes ready-to-use Simple2™ Enzymatic Detergent and an endoscopic cleaning pad in a stand-up pouch.



SKU	OPTIONS	SIZE	UNITS/CASE
IMPORTANT: The following kits are intended for use on small diameter scopes in SPD only.			
EP-2	Endoscopic Cleaning Pad	70 ml	100 Kits
EP-3	Endoscopic Cleaning Pad	100 ml	100 Kits



SKU	OPTIONS	SIZE	UNITS/CASE
EP-4	Endoscopic Cleaning Pad	200 ml	70 Kits
EP-4B	EP4 Kit with Button Brush	200 ml	70 Kits
EP-4G	Gemini Pads	200 ml	70 Kits



SKU	OPTIONS	SIZE	UNITS/CASE
EP-6	Endoscopic Cleaning Pad	500 ml	30 Kits

The First Step Add Water Kit

Includes concentrated Simple2™ Enzymatic Detergent and an endoscopic cleaning pad in a resealable stand-up pouch.

Add water to the Kit when you are ready to use it. The concentrated detergent can be diluted in the pouch, no additional mixing container is required. The Add Water Kit maximizes storage space and reduces freight costs.



Tabletop Kit Holders
 Order#: FS74
 Size: 7"x 4"
 Order#: FS94
 Size: 9"x 4"



Before Water Water Added



SKU	OPTIONS	SIZE	UNITS/CASE
EP-6W	Endoscopic Cleaning Pad	500 ml	120 Kits

Keep Kits Within Reach



Wall Rack
 Order#: WR95M
 Size: 12"H x 7"W x 7"D